


ErgoStar® – catheter mounts and airway connectors

Dräger stands for modern ventilation treatment of high quality and safety when it comes to the care and treatment of patients. ErgoStar® catheter mounts and airway connectors are one part in the line of breathing system equipment from Dräger.



Dräger consumables are tested and manufactured to meet the highest standards in medical technology, and their functionality is convincing. Benefits of ErgoStar are:

- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages
- Extendible product provides a lightweight and flexible solution
- Double swivel connector increases ease of use
- Double cap to maintain PEEP during suctioning or bronchoscopy
- Each catheter mount and airway connector is equipped with a red safety cap so that the system is completely closed and protected until directly prior to use on the patient
- 22F/15M connector on device side as a two-in-one solution
- Standardized connectors for safe connection to other components
- Luer lock port for gas sampling
- Disposable convenience enhances workflow
- PVC-free design of gas leading components eliminates risk of softeners such as DEHP



D-9378-2009

ErgoStar® family

ORDER INFORMATION

ErgoStar CM 40

Part no.: MP01840 Packing unit (pcs.): 50

ErgoStar CM 50

Part no.: MP01850 Packing unit (pcs.): 50

ErgoStar CM 60

Part no.: MP01860 Packing unit (pcs.): 50

ErgoStar AC 90

Part no.: MP01890 Packing unit (pcs.): 100

ErgoStar CM 45

Part no.: MP01845 Packing unit (pcs.): 50

ErgoStar CM 55

Part no.: MP01855 Packing unit (pcs.): 50

ErgoStar AC 95

Part no.: MP01895 Packing unit (pcs.): 100



MT-1922-2008

ErgoStar® CM 40

SPECIFICATIONS

	ErgoStar CM 40	ErgoStar CM 45	ErgoStar CM 50	ErgoStar CM 55	ErgoStar CM 60	ErgoStar AC 90	ErgoStar AC 95
Hose length	10 cm	15 cm	10–16 cm	15 cm	10 cm	–	–
Internal volume	26 mL	34 mL	25–44 mL	31 mL	16 mL	7 mL	8 mL
Resistance	0.21 mbar at 60 L/min	0.21 mbar at 60 L/min	0.24 mbar at 60 L/min	0.09 mbar at 60 L/min	0.17 mbar at 30 L/min	0.21 mbar at 60 L/min	0.11 mbar at 60 L/min
Leakage	<50 mL/min at 60 mbar	<50 mL/min at 60 mbar	<50 mL/min at 60 mbar	<25 mL/min at 60 mbar	<50 mL/min at 60 mbar	<50 mL/min at 60 mbar	<25 mL/min at 60 mbar
Compliance	<1 mL/mbar at 60mbar	<1 mL/mbar at 60 mbar	<1 mL/mbar at 60 mbar	<1 mL/mbar at 60 mbar	<1 mL/mbar at 60 mbar	<1 mL/mbar at 60 mbar	<1 mL/mbar at 60 mbar
ISO connectors	22M/15F – 22F/15M	22M/15F – 22F/15M	22M/15F – 22F/15M	22M/15F – 22F/15M	22M/15F – 15M	22M/15F – 15M	22M/15F – 15M
Product	PVC-free Latex-free	PVC-free Latex-free	PVC-free Latex-free	PVC-free Latex-free	PVC-free Latex-free	PVC-free Latex-free	PVC-free Latex-free

Clean-room manufactured according EN ISO 14644-1:1999



MT-1924-2008

ErgoStar® CM 45



MT-1919-2008

ErgoStar® CM 50



MT-1925-2008

ErgoStar® CM 55



D-9372-2009

ErgoStar® CM 60



MT-1927-2008

ErgoStar® AC 90



MT-1926-2008

ErgoStar® AC 95

HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA

Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 436 24 762
Fax +971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Dräger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:

Dräger Medical GmbH
23542 Lübeck, Germany
The quality management system at
Dräger Medical GmbH is certified
according to ISO 13485, ISO 9001
and Annex II.3 of Directive 93/42/EEC
(Medical devices).