

Dräger SpO₂ Sensors

Reusable and Disposable SpO₂ Sensors.



- Dräger SpO₂ Sensors are compatible with Dräger Infinity® Patient Monitors.
- Dräger SpO₂ Sensors are validated for all Dräger Infinity Patient Monitors.
- Dräger SpO₂ Sensors are compatible with the Nellcor® Oxisensor II.



Dräger SpO₂ Sensors

TECHNICAL DATA

LED specification red 660 nm typical 1.8 mW typical

Infrared 905 nm typical 2.0 mW typical

SpO₂ Accuracy adults or pediatrics: 70-100% ± 2 digits

SpO₂ Accuracy in neonates: 70-100% ± 3 digits

The accuracy of the pulse rate depends on the monitor.

ORDER INFORMATION

Dräger Sensors

	Order number
SpO ₂ Sensor, Adults, Reusable	MS 13 235
SpO ₂ Sensor, Adults, Vinyl, Disposable	MS 16 449
SpO ₂ Sensor, Adults, Foam, Disposable	MS 16 445
SpO ₂ Sensor, Pediatrics, Vinyl, Disposable	MS 16 448
SpO ₂ Sensor, Pediatrics, Foam, Disposable	MS 16 444
SpO ₂ Sensor, Infants, Fabric, Disposable	MS 16 447
SpO ₂ Sensor, Neonates, Fabric, Disposable	MS 16 446

Adapter cable

To connect Dräger/ Nellcor SpO ₂ Sensors to MultiMed 5/6, MultiMed 12 Pod and NeoMed Pod	
Extension cable SpO ₂ cable, 1 m	33 68 433
Extension cable SpO ₂ cable, 2 m	33 75 834
To connect Dräger/ Nellcor SpO ₂ Sensors to MultiMed <i>plus</i> , MultiMed <i>plus</i> OR Pods and M300	
SpO ₂ cable, 1.2 m	MS 18 683
SpO ₂ cable, 3.0 m	MS 17 330



MT-8237-2005

SpO₂ Sensor, Pediatrics, Foam, Disposable



MT-8242-2005

SpO₂ Sensor, Adults, Foam, Disposable



MT-8238-2005

SpO₂ Sensor, Neonates, Fabric, Disposable

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA

Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel + 971 436 24 762
Fax + 971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Dräger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com